



TheraSwitch—

Your mobile assistant for
dysphagia assessment and
treatment during home visits



Without a safe patient position, examination and therapy become more difficult

During home visits, you often encounter situations that make your work difficult: The patient is sitting too low, on a surface that's too soft, and without adequate support. Or they're reclining in a position that makes it hard for you to reach them.

Before you can start an assessment or therapy exercises, you must first find a position that allows you to work effectively.

An optimal position is especially crucial when treating dysphagia.

During FEES examinations, the head must remain still and stable so that the swallowing process can be reliably assessed.

In therapy you also require a stable and relaxed starting position so that you can perform exercises effectively.

In your practice, you can create the right environment for your treatment. During home visits, you have to work with whatever is available.

The challenge:

Part of your attention is focused not on diagnosis and therapy, but on stabilization and balance.

But your goal remains the same:

You want to perform examinations and treatments in a stable position that provides your patients with secure support and allows you to work with full concentration.

Our solution:

The TheraSwitch—ready wherever care takes place.

“The TheraSwitch provides significantly more flexibility in therapy. Even in difficult and complex settings, it ensures that assessments and therapy exercises can be carried out safely and consistently, and it gives you greater independence when it comes to therapy options.”

Sören Paul, MSc,
Therapy Scientist,
Specialist Therapist for Dysphagia, Certified FEES Instructor (DGN, DSG, DGG)

Practical support for everyday therapy

The TheraSwitch wasn't developed in theory—it was shaped by real-world experience. As part of training sessions at the Dysphagiezentrum—Zentrum für Diagnostik und Therapie bei Schluckstörungen (Center for the Diagnosis and Treatment of Swallowing Disorders)—therapists shared their experience and provided input on requirements and usability throughout the development process. The result is a portable assistant that supports you exactly where things used to be most difficult.

The TheraSwitch gives you something that's often missing during treatments:

Stability. Safety. Confidence.

- > A stable head and torso position
- > Less postural effort, less exertion
- > A reproducible starting position for diagnosis and therapy

As a tested and certified medical device manufactured in Germany, the TheraSwitch is a reliable companion for your daily work routine.

"This portable assistant for assessment and therapy should be an indispensable part of every dysphagia specialist's toolkit."

Mirko Hiller,
MSc in Neurorehabilitation,
Owner & Managing Director,
Das Dysphagiezentrum



How to create a stable position for your patients From a chair to a recliner—in just a few simple steps


The backrest, footrest, and headrest can be quickly adjusted.

- > Height-adjustable legs for an ergonomic working posture
- > Armrests can be attached or removed as needed
- > Three different headrests:
 - Shell-shaped cushion for secure support during examinations
 - Flat cushion for relaxed therapy
 - Semi-circular cushion (optional)

Shell-shaped cushion



Flat cushion



Semi-circular cushion



Hygienic & durable

- > Wipeable, disinfectable surfaces
- > Supports up to 120 kg—for a wide range of patients

Mobile & well-designed

- > Weighs only 16 kg—lighter than some wheeled suitcases
- > Large, sturdy wheels—for uneven surfaces and stairs
- > Foldable—fits in standard trunk spaces
- > No electrical or hydraulic components—low-maintenance
- > Convenient accessory bag—everything within easy reach in one place

What changes for you

With the TheraSwitch, it's not just your patients' position that changes—it's the way you work during therapy.

For you as a therapist

- > You'll work more calmly and ergonomically
- > You'll need to hold and balance less
- > You'll be able to focus more fully on your patients and your professional expertise
- > You go home in the evening feeling:
"Today, I was able to focus on my patients and their needs."

For your patients

- > More support and security
- > Less stress thanks to stable positioning
- > Better conditions for assessment and therapy

For your daily practice

Whether it's dysphagia therapy, consistency tests, posture adjustments, or—where available—FEES:

The TheraSwitch provides the foundation so that your work can proceed safely, systematically, and effectively anywhere.



Who the TheraSwitch is designed for

The TheraSwitch is designed for professionals who need safety, support and consistent positioning during home visits and mobile care:

- > **Speech and language therapists**
- > **Therapy practices looking for foldable and compact support**
- > **Dysphagia specialists**
- > **Outpatient teams in neurorehabilitation**
- > **Care facilities and group homes**
- > **Doctors providing mobile care**

Get started in a way that works best for you.

Choose the option that fits your needs: purchase, 0% financing for up to 12 months, flexible payment plans, or trial and rental options.

**Always ready to go:
Fold it up, load it in,
and head to your
next patient.**



The next step: See it, test it, experience it for yourself

Get to know the TheraSwitch.

You have several options:

- > Get a demonstration of the TheraSwitch—online or on-site—
with no obligation.
- > Test it in your own daily work routine and see for yourself what changes.
- > Talk to us about your situation, we'd be happy to help you.

Get a first impression in this short video:

1 minute. No words. Just see how it works.



We're here to help

If you have questions, would like more information, or are considering how our mobile assistant can support you, we're here for you. So you can make an informed decision.

Ergo-Tec GmbH

Peter Prowinsky—Medical Device Consultant

+49 9104 8296997 • sales@ergo-tec.com • www.et-medical.eu